

## Resolute® FM Chromatography Columns



### *Advanced column design enabling true linear scale up*

Resolute FM columns are ideal for chromatographic purification from process development through to manufacturing scale with a wide range of commonly used media. These glass columns, for development and manufacturing applications, are robust and incorporate new features to improve the outcome of column packing events. The industry accepted constant compression seal design used in the successful Resolute DM range is also incorporated in the FM series. The columns are ideal for use with Pall's PKP chromatographic systems providing an ideally suited platform for pilot scale and small volume manufacturing purification operations with a wide range of chromatographic media. Accessory column valves and tubing kits are available to provide a complete solution for your next purification challenge.

#### **Key Features and Benefits**

- ▶ Unique constant compression adjuster seal simplifies column operation
- ▶ Low friction adjuster mechanism enables packing compression to x1.2cf
- ▶ Easy air removal prior to packing for first time packing event
- ▶ Improved safety with an adjuster that enables packing without removal from column
- ▶ Low pressure drop flow path supports high linear velocity up to 700 cm/hr for maximum productivity
- ▶ Simple maintenance of seals and nets with quick and easy removal of both distributor plates without use of special tools
- ▶ Hygienic design suitable for effective sanitization
- ▶ Precision borosilicate glass tube and stainless steel flow path provides excellent chemical compatibility
- ▶ Scale up with 3 column sizes 100, 140 and 200 mm diameter with capacities from 0.78 L to 8.4 L

## Applications

The superior design of Pall's Resolute columns is optimized for purification in both process optimization and small scale manufacturing operations. These columns have proven performance with a wide range of chromatography media and chromatographic modes including: ion exchange, mixed mode, affinity, size exclusion and hydrophobic interaction. Column bed heights up to 35 cm may be prepared with a compression factor of 1.2. For bed heights above 35 cm a longer tube kit is available to support maximum bed heights up to 41 cm.



## Modular Design

Resolute FM columns can be packed and unpacked with the column adjuster remaining supported above the top of the tube. The central adjuster mechanism enables the top distribution cell to travel within the tube and to be raised clear of the tube for packing and unpacking operations. Thus eliminating the requirement to remove and lay down the adjuster and associated tubing, thus improving operator safety and reducing preparation time where connected hoses may be left in place.

The adjuster seal positioned close to the edge of the distribution plate is under constant compression within the tube. There is no need to adjust and fine tune seal compression as a seal is assured every time without operator involvement. This feature avoids problems associated with seals being either under or over compressed. This avoids leakage of fluids past the seal during the lowering of the adjuster during axial bed compression.

Unlike other column designs, the adjuster is always maintained in a relative vertical plane to the column tube by the adjuster plate and does not need to be inserted by

the operator at an angle into the tube to remove air while the seal is formed against the tube wall. Resolute FM columns have a unique feature on the lower frame which allows the column assembly to be tilted from vertical while the adjuster is lowered. This feature has simplified the packing preparation process and results in a reproducible outcome which can be defined in SOP's and avoids awkward handling of heavy column components and attached tubing. This procedure is ideal for concise SOP's and produces reproducible column priming outcomes due the security of the mechanical steps involved in place of manual manipulation of separate column parts.

## Fast Packing and Unpacking

The design of the central adjuster mechanism provides adequate clearance for entry of the slurry directly into the column tube. With a measured volume of slurry in the column the liquid level is then topped up with buffer to a specific bed height (marked on the tube unit). The column is then set in the priming position and the adjuster lowered into the tube by rotating the centrally located handles. As the seals enter the tube the air is displaced and once the seal is made the column is again set vertical before applying the required flow to pack the bed. Once the bed is compressed to a target value the flow is stopped and the adjuster lowered to meet the packed bed. This may be repeated several times dependant on media type.

Unpacking can be performed through the bottom or top of the tube. Typically, the adjuster is raised by rotating the centrally located handles until it is clear of the tube. After rinsing the top distribution plate and net to recover residual media the column contents are gently agitated with buffer to form a slurry before being pumped to a container.




## Performance


Figure 1

Acetone 2% (w/v), injection volume x% CV,  
bed height 20 cm.


FM100



FM140



FM200



## Ease of Use

Resolute FM columns are suitable for sanitization as an empty column (or in a packed condition), as well as having easy access to all seals and supports for fast maintenance requirements. With a minimal number of parts, access for wipe down and the electropolished stainless steel surface finish minimizes entrapment of contaminants and reduces the risk of corrosion. Validation of all column sizes demonstrated that the columns can be effectively sanitized.

## Specifications

Product	Diameter (mm)	CSA (cm <sup>2</sup> )	Adjustable Height Range (mm)	Adjustable Capacity Packed Volume (L)*	Maximum Allowable Pressure (barg/psig)	Net Weight Kg (lb)	Column Height at Minimum Bed (mm)	Width (mm)	Depth (mm)	Height (mm)
Resolute FM	100	79	100 - 460	0.78 - 2.1	5.0/72.5	30/66	1105	385	385	1625
Resolute FM	140	154	100 - 460	1.54 - 4.1	5.0/72.5	45/99	1135	415	415	1655
Resolute FM	200	314	100 - 460	3.14 - 8.4	4.6/66.7	70/154	1135	455	455	1655


\*Based on 70% slurry and compression factor x 1.2

### Temperature

2-35 °C

### Column Termination Sizes

	T1 mm	T2 mm
Diameter (mm) 100	4	25
Diameter (mm) 140	6	25
Diameter (mm) 200	6	25



### Component Weights Kg (lb)

	FM100	FM140	FM200
Column	30 (66.1)	46 (101.4)	62 (136.7)
Adjuster Assembly Complete	8.6 (19.0)	18.57 (41.0)	25.1 (55.3)
Adjuster Distribution Plate Only	1.1 (2.4)	2.2 (4.8)	4.6 (10.1)
Fixed Distribution Plate Only	4.0 (8.8)	6.7 (14.8)	11.9 (26.2)

### Materials of Construction

#### Process Wetted Components

Glass, Borosilicate 3.3; Stainless steel, grade 316 L (EN 1.4404)

Nets: 20 µm – Polypropylene 10 µm – Polyamide

Net support: Polypropylene

#### Seals

EPDM, black, peroxide cured

#### Stand

Stainless steel, grade 316 L (EN 1.4404)

#### Castors

Stainless steel, grade 304

### External Components

#### Column Tube

Glass, Borosilicate 3.3

#### Pressure Envelope Components

Stainless steel, grade 316 L (EN 1.4404)

#### Stand

Stainless steel, grade 316 L (EN 1.4404)

#### Castors

Stainless steel, grade 304

#### Wheels

Acetal

Materials used in the manufacture of the column product contacting components either conform to relevant sections of the FDA Code of Federal Regulations Volume 21, parts 170 to 199, or are non-FDA recognized materials that pass USP Class VI tests for toxicity. All Wetted Elastomers Conform to USP VI or CFR 177.2600

### Applicable Codes and Standards

Design Standard: Steel components/Glass Tube	PD5500 (2009)/BSEN1595 (1997)
Mechanical Design	PED 97/23/EC
P.E.D Category	SEP

## Chemical Combatibility

Substance	Concentration	Stainless Steel	PP	EPDM	Glass	Column overall status	Nylon
Acetic acid	≤ 25%	✓	✓	✓	✓	✓	X
Acetone	50%	✓	✓	✓	✓	✓	✓
Acetonitrile	100%	✓	Softens after ca. 7 days const. contact; effect is reversible	✓	✓	~	✓
Ammonium sulphate	2M	Sl. corrosive to 304 stainless steel. Non-corrosive to 316.	✓	✓	✓	~	✓
Benzyl alcohol	2%	✓	✓	✓	✓	✓	X
Ethanol	100%	✓	✓	✓	✓	✓	✓
Ethanol	40%	✓	✓	✓	✓	✓	✓
Ethanol	20%	✓	✓	✓	✓	✓	✓
Ethyl acetate	100%	✓	✓	✓	✓	✓	✓
Ethylene glycol	50%	✓	✓	✓	✓	✓	✓
Formaldehyde	40%	✓	✓	✓	✓	✓	~
Formic acid	100%	✓	✓	✓	✓	✓	X
Guanidine hydrochloride	8M	X	✓	✓	✓	X	
Hydrochloric acid	20%	X	✓	✓	✓	X	X
Hydrogen peroxide	10%	✓	✓	✓	✓	✓	X
Isopropyl alcohol (IPA, Propan-2-ol)	100%	✓	✓	✓	✓	✓	X
Nitric acid	5%	✓	✓	✓	✓	✓	X
Nitric acid	20%	✓	✓	✓	✓	✓	X
Phosphoric acid	40%	Not recommended	✓	Minor effect	✓	✓	X
Sodium acetate	All conc.'s	✓	✓	✓	✓	✓	~
Sodium bicarbonate	All conc.'s	✓	✓	✓	✓	✓	✓
Sodium chloride	0.5M	At pH < 7 pitting occurs also on storage at any pH > 30 days.	✓	✓	✓	~	✓
Sodium chloride	2M	At pH < 7 pitting occurs also at storage at any pH > 24 hrs	✓	✓	✓	~	✓
Sodium hydroxide	2M	✓	✓	✓	Exposure > 3 months can damage glass surface	~	~
Sodium hydroxide	0.02M	✓	✓	✓	✓	✓	✓
Sodium hypochlorite (200 ppm Cl)	20%	Some discolouration can occur if contact time is > 1 hr	Minor effect	Minor effect	✓	~	X
Sodium nitrate	All conc.'s	✓	✓	✓	✓	✓	✓
Sulphuric acid	≤ 10%	Some minor effect	✓	✓	✓	✓	X
(Tri-(n-butyl) phosphate with 1% Tween)	≤ 3%	✓	✓	✓	✓	~	
Trifluoroacetic acid	≤ 0.1%	Higher conc.'s can cause minor damage to steel	✓	✓	✓	~	
Triton X-100	3%+	✓	✓	✓	✓	✓	
Urea	8M	✓	✓	✓	✓	✓	✓

✓ Compatible

~ Limited compatibility

X Not compatible

## Ordering Information

### Resolute FM Chromatography Column

Item Description	Diameter (mm)	CSA (cm <sup>2</sup> )	Adjustable Height (mm)	Adjustable Capacity Litres	Pressure		Adjustable Capacity Columns	
					psi	bar	Resolute FM	Documentation Pack
FM 100/100-460	100	79	100 - 460	0.78 - 2.1	72	5	FM100A100A20	SFM100A260A99
FM 140/100-460	140	154	100 - 460	1.54 - 4.1	72	5	FM140A100A20	SFM140A260A99
FM 200/100-460	200	314	100 - 460	3.14 - 8.4	69	4.6	FM200A100A20	SFM200A260A99
<b>Tube Kit</b>								
100/350-710	100	79	350-710	2.74 - 3.2	71	5	SFM100A400A99	N/A
140/350-710	140	154	350-710	5.38 - 6.4	72	5	SFM140A400A99	N/A
200/350-710	200	314	350-710	10.99 - 13.0	69	4.6	SFM200A400A99	N/A

NOTES: Column ships with 20 µm nets fitted.

## Spares Kits and Accessories

Spare Kits						
Column Ø (mm)	Fixed & Adjuster Net (10 µm)	Fixed & Adjuster Net (20 µm)	Adjuster Net Support	Fixed Net Support <sup>1</sup>	Wetted Seals <sup>2</sup>	Non Wetted Seals <sup>3</sup>
100	SFM100A150A99	SFM100A160A99	SFM100A350A99	SFM100A360A99	SFM100A300A99	SFM100A340A99
140	SFM140A150A99	SFM140A160A99	SFM140A350A99	SFM140A360A99	SFM140A300A99	SFM140A340A99
200	SFM200A150A99	SFM200A160A99	SFM200A350A99	SFM200A360A99	SFM200A300A99	SFM200A340A99

## Accessories

Accessories					Tools	
Column Ø (mm)	Column Valve Kit <sup>4</sup>	Bottom Valve Mounting Kit <sup>5</sup>	Foot Kit	Tubing Kit <sup>6</sup>	Deep Socket <sup>7</sup>	Stirring Paddle
100	SFM100A270A99	SFM100A320A99	SFM100A310A99	SFM100A420A99	SFM100A330A99	SFM100A390A99
140	SFM140A270A99	SFM140A320A99	SFM100A310A99	SFM140A420A99	SFM140A330A99	SFM100A390A99
200	SFM140A270A99	SFM140A320A99	SFM100A310A99	SFM140A420A99	SFM140A330A99	SFM100A390A99

1. Net support (1) to fit between net and distributor
2. Adjuster seal (1), tube seal (2), inlet tube seal (1)
3. Adjuster nut flange seal (1), inlet spool top seals (2)
4. Two way valve (1), clamp (1), gasket (1) [ 2 valve kits required per column for top and bottom]
5. Bottom elbow (1), clamp and gasket (1), bracket (1) and bolts/washers (2)
6. Reinforced tubing (2 x 150cm) with hose barbs(4), gaskets (4) and clamps (4)
7. Required when dismantling tube assembly

## Documentation

- ▶ Each column is supplied with an Operating and Maintenance Guide
- ▶ CD-ROM with Operating and Maintenance Guide and language versions

## Installation and Validation Documentation

Resolute FM Columns can be ordered with an optional\*\* comprehensive validation and support package. To support validation efforts and efficient commissioning this documentation includes all the pre-shipment test results.

- ▶ Column documentation includes:
  - Setup and test report
  - Materials and compatibility data
  - Comprehensive parts list

\*\*Must be ordered with column at time of purchase

## Chromatography Process Solutions from Pall

Pall offers a comprehensive, versatile and environmentally friendly range of Chromatography Sorbents that greatly simplifies protein purification challenges from lab to process scale. Pall chromatography sorbents can improve the downstream purification accuracy, speed, safety, reliability, and decrease purification cost. Sorbent chemistries include ion exchange, mixed mode, hydroxyapatite, affinity, size exclusion, solvent-detergent removal and HIC.

The PK solution platform and PKP Chromatography System range incorporates proven component technology, delivering accurate precision gradients as well as in-line buffer dilution. This occurs while maintaining the performance necessary to make your chromatography process consistent from batch-to-batch and optimal at pilot and manufacturing scales. Covering a range from 1 < 4000 L/h, Pall's PK and PKP systems ensure robust operation and flexibility for all biochromatography applications. These systems can be used with Resolute chromatography columns as well as Mustang® chromatography membrane adsorber capsules.



Pall LRC chromatography columns have been designed to meet the requirements of most laboratory applications (ion exchange, affinity, size exclusion, mixed-mode chromatography). Four internal diameters, from 10 to 50 mm, and four lengths provide a range of column volumes up to 900 mL and bed heights up to 750 mm. Each column is supplied, assembled, and ready-to-use with tubing and fittings for connection to a standard system. The column body is made of a borosilicate glass tube. The columns are equipped with one adjustable and one fixed plunger, and a screw-lock system, to allow a rapid adjustment of the piston.

The Resolute DM columns are equipped with Pall's proprietary high resolution mobile phase flow path, which consists of a mobile phase flow path within the nozzle valve body, distribution end cells, and a fully flushed piston seal arrangement. The performance of these essential column components is increasingly important for optimal process efficiency, capacity, and peak symmetry when scaling up to larger capacity columns.



Resolute DM columns maintain identical design parameters throughout both the standard column range (from 280 mm to 1200 mm diameter) and the larger capacity columns (up to 2000 mm diameter). Operation at high linear velocities is free from non-linear column effects and assures true linear scalability of column performance, enabling users to scale up simply without method redevelopment or impact on the reproducibility of the packing process.





Life Sciences

**United States**

800.717.7255 toll free (USA)  
516.484.5400 phone  
516.801.9548 fax  
biopharm@pall.com E-mail

**Europe**

+44 (0)23 9230 3303 phone  
+44 (0)23 9230 2506 fax  
BioPharmUK@europe.pall.com E-mail

*Filtration. Separation. Solution.<sup>SM</sup>*

---

**Visit us on the Web at [www.pall.com/biopharm](http://www.pall.com/biopharm)**


**E-mail us at [biopharm@pall.com](mailto:biopharm@pall.com)**

---

**International Offices**

Pall Corporation has offices and plants throughout the world in locations such as: Argentina, Australia, Austria, Belgium, Brazil, Canada, China, France, Germany, India, Indonesia, Ireland, Italy, Japan, Korea, Malaysia, Mexico, the Netherlands, New Zealand, Norway, Poland, Puerto Rico, Russia, Singapore, South Africa, Spain, Sweden, Switzerland, Taiwan, Thailand, the United Kingdom, the United States, and Venezuela. Distributors in all major industrial areas of the world.

The information provided in this literature was reviewed for accuracy at the time of publication. Product data may be subject to change without notice. For current information consult your local Pall distributor or contact Pall directly.

© 2011, Pall Corporation. Pall,  and Mustang are trademarks of Pall Corporation. ® indicates a trademark registered in the USA. *Filtration. Separation. Solution.* is a service mark of Pall Corporation.

2/11, PDF, GN10.4274

USD2776